

## 1. Monitoring

### ► 1. Pulseoximeter PALMCARE



This small unit with a large LCD display provides :

- SpO2
- Adjustable alarm
- Waveform on display
- HR
- Vertical bar graph
- Rechargeable battery with adaptor

### ► 2. Pulseoximeter BPM 200



Main unit with LCD screen

- 1 \* Reusable finger probe
- 1 \* power cable
- 1 \* built in battery
- 1 \* operation manual

### ► 3. Sentry series



- |                 |                       |
|-----------------|-----------------------|
| 04 - SENTRY-N   | NIBP                  |
| 04 - SENTRY-NS  | NIBP-SpO2             |
| 04 - SENTRY-NST | NIBP-SpO2-TRENDSCREEN |

#### **4 . Patient monitor BM3Plus (SpO2 – NIBP – ECG – Resp. Temp.)**



##### **Superior advantages**

- 8.4± color TFT screen with high resolution.
- VGA output connecting to an external monitor directly.
- Visible alarm lamp.

##### **Smart vital technology**

- Arrhythmia analysis : VTAC, VFIB, ASYSTOLIC.
- Pacemaker detection indicated on the ECG waveform display.
- Usable from adult to neonatal patient.

##### **Easy-to-use**

- Trim knob for easy setup and prompt selection.
- Direct keys for frequently used functions.
- Internet upgrade for the latest software installation.

##### **Flexible to transport**

- Battery operation makes the device wireless for 2 hours.
- Direct DC (18VDC) input during vehicle transport.
- Easy mounting with hook-up on hand grip.

##### **Lightweight, Compact, Durable**

- Less than 4kgs including battery.
- Easily fits on a standard hospital bed-rail.
- Withstands the rigors of hospital use and patient transport.

##### **Data management and reporting**

- 24 full hours trend and alarm information are stored.
- Built-in recorder generates real time 3-trace data.
- Central monitoring system up to 16 beds. (Option)

29-72030000	BM3-PLUS MONITOR with build-in printer
29-72080110	BM3-PLUS EXTENSION CABLE
29-72080120	ECG-3CABLE for BM3-PLUS
29-72080140	ECG-3LEADS WIRE for BM3-PLUS
29-72080170	PAPER-400 for BM3-plus
29-72080255	MOBILE CART B-400