

16. HEINE



Heine_laryngoscopen

06-0022100	Macintosh Blade Size 0 FO
06-0022101	Macintosh Blade Size 1 FO
06-0022102	Macintosh Blade Size 2 FO
06-0022103	Macintosh Blade Size 3 FO
06-0022104	Macintosh Blade Size 4 FO
06-0022105	Macintosh Blade Size 5 FO
06-0122860	Handle Fiber Optic
Miller also available	



Heine_stroombronnen



Heine_reservelampjes



Heine_laryngoscoop_handvaten



Heine_beta200_oogspiegel



Heine_beta200_otoscope



Heine_Delta20_dermatoscope



Heine_3S_LED_voorhoofdlamp



Heine_beta_skiascoop



Heine_EN100_basis_uitbouw_tipdispenser_2006



Heine_EN100_basis2



Heine_gamma_E60



Heine_gamma_g5



Heine_gamma_manchetten_1T



Heine_gamma_manchetten_2T



Heine_gamma_stethoscoop



Heine_gamme_XXL



Heine_gamme_XXL



Heine_koudlichtprojector



Heine_MD1000



Heine_loepbrillen



Heine_spatelhouder



Heine_sets_a141



Handle FO

HEINE GAMMA G®-Series
Shockproof and Latex-free Sphygmomanometers




06-M0009235	HEINE G5	STANDARD SPYGHMOMANOMETER PLASTIC
06-M0009237	HEINE G7	STANDARD SPYGHMOMANOMETER ALUMINIUM
06-M0009240	HEINE GP	STANDARD SPYGHMOMANOMETER WITH PUSH-BUTTON VALVE
06-M0009271	HEINE GST	STANDARD SPYGHMOMANOMETER WITH INTEGRATED STETHOSCOPE
06-M0009160	HEINE E60	AUTOMATIC SPYGHMOMANOMETER
06-M0009310	HEINE GAMMA XXL-T	TABLE MODEL DE LUXE
06-M0009311	HEINE GAMMA XXL-W	WANDMODEL LUXE
06-M0009312	HEINE GAMMA XXL-R	RAILMODEL
06-M0009313	HEINE GAMMA XXL-S	STATIEFMODEL
06-M0009307	HEINE	RAILCLAMP FOR XXL-R

Latex-free blood-pressure cuffs


Cuffs for GAMMA G7/G5/GP/GST/E60/XXL LF sphygmomanometers.



- ∴ Latex-free.
- ∴ Fit the most widely-used sphygmomanometers.
- ∴ High-quality velcro cuffs made of wear-resistant material. Physiologically tested.
- ∴ Double-sewn.
- ∴ Washable.
- ∴ Index mark for exact positioning.
- ∴ Clear indication of size and control mark for arm circumference.
- ∴ Extra-thick, highly-flexible tube. For long life.
- ∴ Proven reliability. Absolutely airtight system.




GAMMA G7/G5/GP	Single-tube cuffs (cpl.)			Single-tube bladders		Outer cuffs	
	Size	Circumference/cm		Size		Size	
Neo-nate	6x28	7-10	M-00.09.652	5x15	M-00.09.661	6x28	M-00.09.051
Infant	8x34	8-13	M-00.09.639	6x20	M-00.09.662	8x34	M-00.09.052
Child	10x34	14-21	M-00.09.637	8x20	M-00.09.663	10x34	M-00.09.053
Adult	14x53	22-32	M-00.09.617	11x23	M-00.09.664	14x53	M-00.09.054
Oversize	15x64	33-41	M-00.09.613	13.5x40	M-00.09.666	15x64	M-00.09.056
Thigh	20x86	> 41	M-00.09.615	18x52	M-00.09.667	20x86	M-00.09.057



GAMMA XXL LF	Twin-tube cuffs* (cpl.)			Twin-tube bladders		Outer cuffs	
	Size	Circumference/cm		Size		Size	
Neo-nate	6x28	7-10	M-00.09.653	5x15	M-00.09.671	6x28	M-00.09.051
Infant	8x34	8-13	M-00.09.643	6x20	M-00.09.672	8x34	M-00.09.052
Child	10x34	14-21	M-00.09.641	8x20	M-00.09.673	10x34	M-00.09.053
Adult	14x53	22-32	M-00.09.627	11x23	M-00.09.674	14x53	M-00.09.054
Oversize	15x64	33-41	M-00.09.623	13.5x40	M-00.09.676	15x64	M-00.09.056
Thigh	20x86	> 41	M-00.09.625	18x52	M-00.09.677	20x86	M-00.09.057

*includes male/female connectors. Compatible with similar bosco units.



To fit other manufacturers	Double-tube cuffs (cpl.)			Double-tube bladders		Outer cuffs	
	Size	Circumference/cm		Size		Size	
Neo-nate	6x28	7-10	M-00.09.681	5x15	M-00.09.691	6x28	M-00.09.051
Infant	8x34	8-13	M-00.09.682	6x20	M-00.09.692	8x34	M-00.09.052
Child	10x34	14-21	M-00.09.683	8x20	M-00.09.693	10x34	M-00.09.053
Adult	14x53	22-32	M-00.09.684	11x23	M-00.09.694	14x53	M-00.09.054
Oversize	15x64	33-41	M-00.09.686	13.5x40	M-00.09.696	15x64	M-00.09.056
Thigh	20x86	> 41	M-00.09.687	18x52	M-00.09.697	20x86	M-00.09.057