

08. New – Hospicode Label system

Reference	Description
14 – ML	Standard label – 1 box = 400 labels
14 – MLDESK	Blue outer container for 10 boxes of 400 labels
14 – MLWALL	Wall container for 10 boxes of 400 labels
www.hospicode.com	

Induction Agents	Propofol mg/ml.	Yellow
Tranquillizers	Midazolam mg/ml.	Orange
Tranquillizer Antagonists	Flumazenil micrograms/ml.	Orange
Non-Depolarizing Muscle Relaxants	Atracurium mg/ml.	Red
Depolarizing Muscle Relaxants	Suxamethonium mg/ml.	Red
Muscle Relaxant Reversal Agents	Neostigmine mg/ml.	Red
Opiates	Fentanyl micrograms/ml.	Blue
Opiate Antagonists	Naloxone micrograms/ml.	Blue
Local Anaesthetics	Lidocaine mg/ml.	Grey
Vasopressors	Dobutamine mg/ml.	Violet
SLS101 & SLS102 only	Adrenaline micrograms/ml.	Violet
Hypotensive Agent	Glyceryl Trinitrate mg/ml.	Violet
Anti-Cholinergic Agents	Atropine mg/ml.	Green
Anti-Emetics	Ondansetron mg/ml.	Salmon
Miscellaneous Drugs	Actrapid Insulin units/ml.	White
Miscellaneous Antagonist Agents	Protamine mg/ml.	White