













9. SpO₂ – sensors and adaptor cables

► Compatible Reusable SpO₂ sensors




























































 <p>A520 Series Reusable Finger Clip</p>	
 <p>A120 Series Reusable Small Finger Clip</p>	
 <p>A220 Series Reusable Ear Clip Sensor</p>	
 <p>A320 Series Reusable 'Y' Sensor</p>	
 <p>A900 Series Reusable SoftTip™</p>	
 <p>A800 Series Reusable SoftTip™</p>	































► **Compatible Reusable SpO₂ sensors (example for A520 finger clip sensor)**

Manufacturer	Model	Cable length(m)	Adult Finger sensor
BCI	All	1.2	PR- A520 -1001N
CSI (Criticare Systems intern.)	501, 503, 504, 507, 508, 581, 1100, Poet (602-3A), Scholar	1.2	PR- A520 -1002N
	501, 503, 504, 507, 508, 581, 1100, Poet (602-3A), Scholar	3.6	PR- A520 -1002
Datascope	Accusat, Accutorr 3/4 SAT, Accutorr Plus (Datascope SpO ₂), Passport (Datascope SpO ₂)	1.2	PR- A520 -1004N
	Accusat, Accutorr 3/4 SAT, Accutorr Plus (Datascope SpO ₂), Passport (Datascope SpO ₂)	3.6	PR- A520 -1004
	Expert, Passport (Nellcor® SpO ₂)	1.2	PR- A520 -1011N
	Accutorr Plus (Nellcor® SpO ₂)	1.2	PR- A520 -1011N
Datex-Ohmeda	Datex with adaptor cable	1.2	PR- A520 -1005N
	Datex without adaptor cable	3.6	PR- A520 -1005
	Ohmeda without adaptor cable	3.6	PR- A520 -1017
	Ohmeda with adaptor cable	1.2	PR- A520 -1017N
GE Medical systems	7000 Series, TRAM A.R., TRAM Modular x00SL Series	3.6	PR- A520 -1017
	Eagle, Dash, SOLAR, TRAM Modular x50SL Series	1.2	PR- A520 -1011N
	Hellige Eagle 1000, 3000, SMU, SMK, SMV	1.2	PR- A520 -1011N
	Corometrics 188, 129	1.2	PR- A520 -1011N
	Corometrics 118, 120, 510, 556	1.2	PR- A520 -1011N
	Corometrics 511	3.6	PR- A520 -1017
	Critikon Dinamap Oxytrak, 8280	1.2	PR- A520 -1003N
	Critikon Dinamap 8700 / 8710 / 8720 Plus	3.6	PR- A520 -1003P
	Critikon Dinamap 8700 / 8710 / 8720 Plus, Oxyshuttle	1.2	PR- A520 -1043N
	Critikon Dinamap 9700 / 9710 / 9720 Plus, Compact, MPS, Pro 300/400, Select	1.2	PR- A520 -1011N
Nellcor	All with adaptor cable	1.2	PR- A520 -1011N
	All without adaptor cable	3.1	PR- A520 -1011
Nihon Kohden	Nihon Kohden without adaptor cable	3.0	PR- A520 -1013
	Nihon Kohden with adaptor cable	1.2	PR- A520 -1013N
	AL-900P Module; BSM-1102, 9510; OGS-2000	1.2	PR- A520 -1011N
	BSM-4104A, 4114A	1.2	PR- A520 -1011N
Nonin	8500, 8600, 8700, 8800	1.2	PR- A520 -1014N
	8604	3.0	PR- A520 -1014L
Novamatrix	510, 511, 505, 515B/C, 520, 7000, 7100, Oxypleth	3.0	PR- A520 -1016
	510, 511, 505, 515B/C, 520, 7000, 7100, Oxypleth	1.2	PR- A520 -1016N
	500	3.0	PR- A520 -1015
Palco	All	1.2	PR- A520 -1018N
Philips Medical Hewlett Packard	78352C, 78354C, 78834C, M1020A, M1025B, M1205A, M1350B, M1722A/B, M1732A/B, M2475B	3.6	PR- A520 -1006
	78352C, 78354C, 78834C, M1020A, M1025B, M1205A, M1350B, M1722A/B, M1732A/B, M2475B	1.2	PR- A520 -1011N
Spacelabs Medical	90419, 90501, 90651A-03	3.0	PR- A520 -1015
	90419-01, 90651A-07	3.0	PR- A520 -1016
	90351-D/G, 90465, 90466, 90467, 90489, 90651-07/08	3.0	PR- A520 -1016
	90351-D/G, 90465, 90466, 90467, 90489, 90651-07/08	3.0	PR- A520 -1020
	90309 PC Scout	3.0	PR- A520 -1010PC
	90351-0/6, 90465, 90466, 90467, 90489, 90651A-08, 90651A, IM77	1.2	PR- A520 -1011N
	90496 Ultraview (Nellcor® SpO ₂)	1.2	PR- A520 -1011N
	90496 Ultraview (Novamatrix® SpO ₂)	3.0	PR- A520 -1016

For A120, A220, A320, A800, A900 adjust reference code

► **Adaptor cables**

Part No.	Cable Length	Sensor Side	Monitor Side
CB-A400-1001A BCI	3.0m	100 	184  
CB-A400-1001E BCI	2.4m	100 	101  
CB-A400-1002A Criticare	3.0m	100 	177  
CB-A400-1002E Criticare	3.0m	100 	105  
CB-A400-1004A Datascope	3.0m	100 	164  
CB-A400-1004B Datascope	3.0m	100 	164  
CB-A400-1005A Datex-Ohmeda	3.0m	100 	011  
CB-A400-1006A HP-Philips	3.0m	100 	025  
CB-A400-1006V HP-Philips	3.0m	088  	025  
CB-A400-1006VM HP-Philips	1.98m	088  	087  
CB-A400-1006VN HP-Philips	3.0m	100 	087  
CB-A400-1011D4 Nellcor	1.2m	100 	110  
CB-A400-1011D8 Nellcor	2.4m	100 	110  
CB-A400-1011D90 Nellcor	2.4m	100 	110  
CB-A400-1011DO Nellcor	3.0m	100 	097  
CB-A400-1011E Nellcor	2.4m	100 	103  
CB-A400-1011E90 Nellcor	2.4m	100 	103  
CB-A400-1011ES Nellcor	1.2m	100 	103  
CB-A400-1011OE Nellcor	3.0m	100 	073  

Part No.	Cable Length	Sensor Side	Monitor Side
CB-A400-1011PA Nellcor	3.0m	100 	180 
CB-A400-1011SC Nellcor	3.0m	100 	073 
CB-A400-1013A Nihon Kohden	3.0m	100 	023 
CB-A400-1016A Novamatrix	2.4m	100 	042 
CB-A400-1016E Novamatrix	3.0m	055 	042 
CB-A400-1017A Ohmeda	2.4m	100 	046 
CB-A400-1020A Spacelabs	3.0m	100 	050 
CB-A400-1020B Spacelabs	3.0m	055 	050 
CB-A400-1024A Siemens	3.0m	100 	031 
CB-A400-1027A GEMS	3.0m	187 	179 
CB-A400-1040A GE	1.8m	100 	052 
CB-A400-1049AN Nellcor-Masimo	3.0m	100 	073 
CB-A400-1102A GEMS	3.0m	100 	064 
CB-A400-1103A Spacelabs	3.0m	100 	074 
CB-A400-1104A Spacelabs	3.0m	055 	074 


► **Original Nellcor sensors**

Durasensor® Finger Clip sensor (Oximax)	
Reference	Description
NU-DS-100A	Pulse oximeter, Nellcor #DS100A, Durasensor, adult



Durasensor® Finger Clip sensor

Oxisensor® II Adhesive Disposable Sensors (Oximax)		
Reference	Description	U/M
NU-S-Dmax-P	Original Nellcor Oximaxsensor, pediatric	Box/24 pcs
NU-S-Dmax-A	Original Nellcor Oximaxsensor, adult	Box/24 pcs
NU-S-Dmax-I	Original Nellcor Oximaxsensor, infant	Box/24 pcs
NU-S-Dmax-N	Original Nellcor Oximaxsensor, neonate	Box/24 pcs



Oxisensor® II Adhesive sensors