

## 11. Breathing systems (Jackson Reese and Watersset)

### ► Adult

<b>Watersset Mapleson C System Latexfree and autoclavable</b>	
Reference	Description
01-0400250LF	System complete latexfree
<b>Parts</b>	
CHM-5003020	Gasbag 2 liter
NU-R-V12200	Gasbag adapter 15F/22M to 22F
01-0400130	Waters expiratory valve with exhaust
NU-R-V12000	Mask elbow 22F to 22M/15F
NU-R-V30002	Feed Mount
23-tubing	Silicone tubing



Watersset Mapleson C System

### ► Child

<b>Jackson Rees system C – Pediatric Latexfree and autoclavable</b>	
Reference	Description
01-0400440LF	System complete latexfree
<b>Parts</b>	
CHM-5003005	Gasbag 0,50 liter
NU-R-V12110	Gasbag adaptor 15F/22M x2
01-0400040	Jackson Rees expiratory valve with exhaust
NU-R-V65002	Mask elbow 22M/15F to 15M
NU-R-V30002	Feed Mount
23-tubing	Silicone tubing



Jackson Rees System C