



Meditron

Catalogoog



info@meditron.be



03/776.22.95

2021-2022

Infrarood thermometer

3



Art. : 04-Therm

Specifications



Measurement Method	Infrared measurement, Non-Contact measurement
Measurement distance	2 to 3cm from the measuring point
Measurement point	The measuring point is the temple area
Measurement range	Body temperature : 22°C to 42.5°C (71.6°F to 108.5°F) Accuracy : ±0.2°C (36°C~39°C), Other Accuracy : ±0.3°C Object temperature : 10°C to 80°C (50°F to 176°F) Accuracy ±2°C
Response time	Less than 2 seconds continuous measurement
Resolution	0.1°C
Memory	Automatically body temperature data (maximum 32) stored
Back-Light	If measurement complete, lighted 3 seconds
Automatic Power off	Auto Power-Off (after 1minute)
Accessory	FS-700 1set, Manual, Battery(AAA type) 2ea
Power supply	DC3V(AAA 1.5V Battery 2ea)



info@meditron.be

03/776.22.95



Meditron

Infrarood thermometer

4



- 1 LCD display
- 2 Start button
- 3 Sound button

- 4 MODE button
- 5 MEM button
- 6 Strap position

- 7 Probe
- 8 Battery cover

Contactloos meten, dus geen overdracht van virussen of bacteriën!

De contactloze koortsthermometer Thermofinder II van Hubdic wordt geplaatst op een afstand van 2 - 3 cm voor de Tempels van de patiënt gehouden en al na 1 - 2 seconden de lichaamstemperatuur verschijnt op het display.

- Meting van de lichaamstemperatuur, zelfs bij slapende patiënten
- Snelle meting binnen 1 - 2 seconden
- Schermresolutie in stappen van 0,1 ° C
- Eenvoudig af te lezen display met achtergrondverlichting
- Ergonomisch ontwerp
- Hoge meetnauwkeurigheid
- Opslag van de laatste 32 meetresultaten
- Batterijbesparende werking door automatische uitschakeling na 60 sec.
- Omschakelbaar voor metingen van de omgevingstemperatuur
- Gebruik via AAA - standaardbatterijen
- Koortsalarm



info@meditron.be

03/776.22.95



Meditron